



Recombinant *Oryza sativa* subsp. japonica B3 domain-containing protein Os06g0107800 (Os06g0107800, LOC_Os06g01860)

Product Code	CSB-YP722785OFG
Abbreviation	Os06g0107800, LOC_Os06g01860
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VS55
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MATIVAWESR NLQLQGGGGG HGGGGGGGGG ERREYMF EKV VTPSDVGKLN RLVVPKH YAE KYFPLGPAAR TSPAGTVLCF EDARGGDSTW RFRYSYWSSS QSYVITKWS RYVRDKRLAA GDTVSF CRAG ARLFIDCRKR AASVSSSSLV PPALIKVQLP PSRPVVDEEE AACGRRCLRL FGVDLQLRAD ASPALDLQL
Source	Yeast
Target Names	Os06g0107800
Protein Names	Recommended name: B3 domain-containing protein Os06g0107800
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.