



Recombinant Human Muscular LMNA-interacting protein (MLIP)

Product Code	CSB-YP722793HU
Abbreviation	MLIP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VWP3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MELEKREKRS LLNKNLEEKL TVSAGGSEAK PLIFTFVPTV RRLPTHTQLA DTSKFLVKIP EESSDKSPET VNRSKSNDYL TLNAGSQQER DQAKLTCPSE VSGTILQERE FEANKLQGMQ QSDLFKAIEYV LIVDSEGEDE AASRKVEQGP PGGIGTAAVR PKSLAISSSL VSDVVRPKTQ GTDLKTSSHP EMLHGMAPQQ KHGQQYKTKS SYKAFAAIPT NTLLEQKAL DEPAKTESVS KDNTLEPPVE LYFPAQLRQQ TEELCATIDK VLQDLSMHS SDSPSRSPKT LLGSDTVKTP TTLPRAAGRE TKYANLSSPS STVSESQLTK PGVIRPVPVK SRILLKKEEE VYEPNPFISKY LEDNSDLFSE QDVTVPPKPV SLHPLYQTKL YPPAKSLLHP QTLSHADCLA PGPFSHLSFS LSDEQENSHT LLSHNACNKL SHPMVAIPEH EALDSKEQ
Source	Yeast
Target Names	MLIP
Protein Names	Recommended name: Muscular LMNA-interacting protein Alternative name(s): Muscular-enriched A-type laminin-interacting protein
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.