



Recombinant *Bacillus clausii* Probable DNA-directed RNA polymerase subunit delta (rpoE)

Product Code	CSB-EP722808BAAC
Abbreviation	rpoE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5WB42
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus clausii</i> (strain KSM-K16)
Purity	≥85% (SDS-PAGE)
Sequence	MDLQSLTKEE LNEMSMVEAA FIVMANKRQP FNFYDLLEEV AKLKGISETE MDERISFLYT DLNLDGRFLT LGDNTWGLRS WYPHEQAEIE ITHSAKPAKR RKSGDDDGEE EYDLVDELEE DEFEDLEDEL DELADEEDQD GERDDYDEKD DALPGFDEED GQTLDEED
Source	<i>E.coli</i>
Target Names	rpoE
Protein Names	Recommended name: Probable DNA-directed RNA polymerase subunit delta Alternative name(s): RNAP delta factor
Expression Region	1-168
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.