



Recombinant *Bacillus clausii* Homoserine kinase (thrB)

Product Code	CSB-YP722815BAAC
Abbreviation	thrB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5WDT6
Product Type	Recombinant Protein
Immunogen Species	Bacillus clausii (strain KSM-K16)
Purity	≥85% (SDS-PAGE)
Sequence	MAFSITVPAS TANLGPGFDS IGLALNRYLH LEVEQASKWS FICSSPGLEN IGTDNLVTKA ATFASKEWGK VLPPCQVVMK NDIPLSKGFG SSAAAI VAGI ELAASVCGQF ASKA EKARLA SLWEGHPDNV AASIYGGLVI GTHSEESTE I LHVEAPNIDL VALIPSKTLA TKKARAVLPD MLSYKEAVQA SSVANVLAAA LVKQDWELVG KMMLADRFHQ PYRQQLIPHL PDVCAYA QAD EDAYGAALSG AGPIILCLVK EGKG EAFASR LHKQFPACRA EVMRPAVQGS AVSIAQA
Source	Yeast
Target Names	thrB
Protein Names	Recommended name: Homoserine kinase Short name= HK Short name= HSK EC= 2.7.1.39
Expression Region	1-297
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.