



Recombinant Legionella pneumophila Putative phosphotransferase lpp1485 (lpp1485)

Product Code	CSB-MP722944LAAB
Abbreviation	lpp1485
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5X535
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila (strain Paris)
Purity	≥85% (SDS-PAGE)
Sequence	MKRYVFMLSD GTGITAETLG NSLITQFENI QFEKITIPYI DSTHRAESVV LRINQCFSEQ GTKPLVFM TL VDPEIRQAIK KAHACVFDLF SIFIGPLENE LEEKSSYTVG RTHGVANVKS YSHRIEADF ALSHDDGIKT RGYDKADIIL IGVSRGKTP SCLYMALQYG ILAANYPFTE EDLVGFRLPE VLRPYKPKLF GLTIDAQRLQ QIRSERRPNS KYASAEQCRL EVTEVEAMYQ RENIPYINST KYSIEEISTK VLAIALGLQRK I
Source	Mammalian cell
Target Names	lpp1485
Protein Names	Recommended name: Putative phosphotransferase lpp1485 EC= 2.7.-.-
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.