



Recombinant Protein RecA (recA)

Product Code	CSB-EP722971SMZ-B
Abbreviation	recA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5X9I0
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MAKKLKKNEE ITKKFGDERR KALDDALKNI EKDFGKGVAM RLGERAEQKV QVMSSGSLAL DIALGAGGYP KGRIIEIYGP ESSGKTTVAL HAVAQAQKEG GIAAFIDAEH ALDPAYAAAL GVNIDELLS QPDSGEGGLE IAGKLIDSGA VDLVVVDVA ALVPRAEIDG DIGDSHVGLQ ARMMSQAMRK LSASINKTKT IAIFINQLRE KVGVMFGNPE TTPGGRALKF YASVRLDVRG TTQIKGTGDQ KDSSIGKETK IKVVKNKVAP PFKVAEVEIM YGEGISRTGE LVKIASDLDI IQKAGAWFSY NGEKIGQGSE NAKRYLAEHP ELFDEIDRVK RVKFGLLLEES EESAMAVAS EETDDLALDL DNGIEIED
Source	E.coli
Target Names	recA
Protein Names	Recommended name: Protein RecA Alternative name(s): Recombinase A
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>