



Recombinant Foldase protein prsA 2 (prsA2)

Product Code	CSB-EP722973SMZ-B
Abbreviation	prsA2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5X9P6
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	CQSSHNNT KLVSMKGDTI TVSDFYNETK NTELAQKAML SLVISRVFET QYANKVSDKE VEKAYKQTAD QYGTSFKTVL AQSGLTPETY KKQIRLTKLV EYAVKEQAKN ETISKKDYRQ AYDAYTPTMT AEIMQFEKEE DAKAALEAVK AEGADFAAIA KEKTIAADKK TTYTFDSGET TLPAEVVRAA SGLKEG NRSE IITALDPAIS KRTHYHIKVT KKATKKADWK AYQKRLKDII VTGKCLKDPDF QNKVIAKALD KANVKIKDKA FANILAQFAK PNQKQPAQK
Source	E.coli
Target Names	prsA2
Protein Names	Recommended name: Foldase protein prsA 2 EC= 5.2.1.8
Expression Region	23-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.