



# Recombinant tRNA modification GTPase MnmE (mnmE)

<b>Product Code</b>	CSB-MP722995SMZ
<b>Abbreviation</b>	mnmE
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5XCB7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSITKEFDTI TAISTPLGEG AIGIVRLSGT DALAIAQSVF KGKNLEQVAS HTINYGHIID PKTGTIIDEV MVSVMLAPKT FTRENVVEIN THGGIAVTNE ILQLLRQGA RMAEPGEFTK RAFLNGRVDL TQAEAVMDII RAKTDKAMTI AVKQLDGSL S QLINDTRQEI LNTLAQVEVN IDYPEYDDVE EMTTALLREK TQEFQSLEN LLRTAKRGKI LREGLSTAI GRPNVGKSSL LNNLLREDKA IVTDIAGTTR DVIEEYVNIK GVPLKLVDTA GIRETDDLVE QIGVERSKKA LQESDLVLLV LNASEKLTQ DRALLNLSQD SNRIILLNKT DLEQKIELEQ LPDDYIPI SV LTNQNINLIE DRINQLFFDN AGLVEQDATY LSARHISLI EKAVQSLEAV NDGLALGMPV DLLQVDLTRT WEILGEITGD AAPDELITQL FSQFCLGK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	mnmE
<b>Protein Names</b>	Recommended name: tRNA modification GTPase MnmE EC= 3.6.-.-
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.