



Recombinant Transcriptional repressor AdcR (adcR)

Product Code	CSB-YP723015SMZ
Abbreviation	adcR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XEA3
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MGTLEKKLDN LVNTILLKAE NQHELLFGAC QSDVKLTNTQ EHILMLLSQQ RLTNTDLAKA LNISQAAVTK AIKSLVKQDM LAGTKDTVDA RVTYFELTEL AKPIASEHTH HHDETLNVYN RLLQKFSAKE LEIVDKFVTV FAEELG
Source	Yeast
Target Names	adcR
Protein Names	Recommended name: Transcriptional repressor AdcR
Expression Region	1-147
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.