



Recombinant Chicken OTU domain-containing protein 6B (OTUD6B)

Product Code	CSB-BP723156CH
Abbreviation	OTUD6B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZIP6
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MEGSEDEEAE AGGPLQQLVK RQRREKRELQ AKIQGMKNAV PKNDKKRRKQ LAEEVAKLEA ELEQKHKEEL KQLKEAMPEQ NKIDSIADGV ANFELEGREQ QIQHPRISKA QKRREKKAAL EKEREERIAE AEIENLTGAR HLESQKLASL LAARHLEIKQ IPSDGHCMYR AIEDQLKDDH NSWTVATLRN QTAKYIHSHF DDFLPFLTNP NTGDMYSKEE FEKYCDDIAN TAAWGGQLEL RALSHILQTP IEVVQMDSPS IIVGEEYSKG PIILVYMRHA YGLGEHYNSV KLLTDATTEN GS
Source	Baculovirus
Target Names	OTUD6B
Protein Names	Recommended name: OTU domain-containing protein 6B
Expression Region	1-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.