



Recombinant Chicken Coiled-coil domain-containing protein 61 (CCDC61)

Product Code	CSB-MP723160CH
Abbreviation	CCDC61
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZJ07
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	MEEPRCLQAN CFFRGKEHSI QLSVSQAVLE VEVEERRSTK RWRGHFDAAS VEDLTRKTGN FKQFGIFCSM LEAALMKSSSE AVSLELLTYG DLEALRCCKV GVATRVPPST SPLSSKRYLI LVYCVDFDRI HYPLPLPYIG EADMAALRRL VQEQQDELAQ LRDELRRRAQQ EVWRLEDERL RDKAWHQEQE RRMTKELTEV KAAEKMLRAH VKTLTAEALV CRKGRSTSAT APGCPQDRRR SNSRDSRSSS QGRLPPRSPS PAGSRPPRFN PTAFVRSREQ RRQEAELRRQ KLPRGTVSSG DNCRRRGRRS SSAESLQSRR SAQSSGSEAD TCPQRRRGLG GPSTRSPLSA SSCNSTSVAS HPNRRGCKQHG KENLGTEPSA TLSEIDARLQ ALQAYISTLG THM
Source	Mammalian cell
Target Names	CCDC61
Protein Names	Recommended name: Coiled-coil domain-containing protein 61
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.