



Recombinant Chicken Endophilin-B2 (SH3GLB2)

Product Code	CSB-MP723163CH
Abbreviation	SH3GLB2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5ZJ81
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MDFNVKKLAS DAGVFFSRAM QFTEEKLQQA EKTELDHFE NLLARADSTK NWTEKILRQT EVLLQPNPSA RVEEFLYEKL DRKVPSRVTN GELLAQYMTE AANDFGPGTP YGKTLIKVGE TQRRLGAAER DFIHSASINF LTPLRNFLEG DWRTISKERR ILQNRRLDLD ACKARLKKAK AAEAKAAAVP DFQETRPRNY VLSASASALW SDEVEKAEHE LRLTQTEFDR QAEVTRLLLE GISSTHVNHL RCLHEFVESQ TNYAQCQYQY MLDLQKQLGR FSGTFVGNAE STSPPPAAAS PPAVAATTLV AVPTIPVVPT IVGVPNTVAE SVLNPNEVKP PASGTRKARV LYDYEAADST ELALLADEMI TVYSLPGMDP DWLIGERGNQ KGKVPVTYLE LLS</p>
Source	Mammalian cell
Target Names	SH3GLB2
Protein Names	Recommended name: Endophilin-B2 Alternative name(s): SH3 domain-containing GRB2-like protein B2
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>