



Recombinant Chicken Muscleblind-like protein 1 (MBNL1)

Product Code	CSB-BP723183CH
Abbreviation	MBNL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZKW9
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MAVSVTPIRD TKWLTLEVCR EFQRGTCSRP DTECKFAHPS KSCQVENGRV IACFDSLKGR CSRENCKYLH PPPHLKTQLE INGRNNLIQQ KNMAMLAQQM QLANAMMPGA PLQPVPMFSV APSLATNASA AFNPYLG PVS PGLVPAEILP TAPMLVAGNP GVPVPA AAAAAA AAQKLMRTDR LEVCREYQRG NCNRGENDCR FAHPADSAMI DTNDNTVTVC MDYIKGRCSR EKCKYFHPPA HLQAKIKAAQ YQVNQAAAAQ AAATAAAMGI PQAVLPPLPK RPALEKTNGA TAVFNTGIFQ YQQALANMQL QQHTAFLPPV PMVHGATPAT VSAATTSATS VPFAATATAN QIPIISAEHL TSHKYVTQM
Source	Baculovirus
Target Names	MBNL1
Protein Names	Recommended name: Muscleblind-like protein 1
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.