



Recombinant Chicken Protein CDV3 homolog (CDV3)

Product Code	CSB-EP723188CH
Abbreviation	CDV3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZLH9
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	AETEERSLD DFFAKRDKKK RKEKSGRSAA AASTASSASS AASSAAGAAA AAGGPRPTDG GGSGAAASSS GSAKSKEEDD WKEFEQKEEI DYSGLRVQSM QISEKEDDES EKREEPSDNW EETSGAVDRS SGPWNKSAPA PAPVVETIVT EPPEPVQTGG VYRPPGAREG GRPRRAQQGP PEIYSDTQFP SLQSTAKLVD SRKY
Source	E.coli
Target Names	CDV3
Protein Names	Recommended name: Protein CDV3 homolog
Expression Region	2-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.