



Recombinant Chicken cAMP-regulated phosphoprotein 19 (ARPP19)

Product Code	CSB-YP723196CH
Abbreviation	ARPP19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZLY8
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	SAESPEPAS AEEQKEMEDK VISPEKAEAA KLKARYPHLG QKPGGSDFLR KRLQKGQKYF DSGDYNMAKA KMKNKQLPTA APDKTEVTGD HIPTPQDLQP RKPSLVASKL AG
Source	Yeast
Target Names	ARPP19
Protein Names	Recommended name: cAMP-regulated phosphoprotein 19 Short name= ARPP-19
Expression Region	2-112
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.