



Recombinant Chicken Rho GTPase-activating protein 15 (ARHGAP15)

Product Code	CSB-EP723208CH-B
Abbreviation	ARHGAP15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZMM3
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MERSTTSDTA SEKPNPSHST GAVQMRIKNA NSHHDRLSQS KSMILSDNVK VLEPINRHRH NHSQHNLTLA DIISTPDHTV VEKEGYLLKA KIADGGKCLR KNWSTSWIVL TSRKMEFYKE SKQPALANLK PGYKPECVDL CGAHIEWTPE KSSRKNVFIQ TTVSGNEFLL QSDIDFLILD WFHAIKNAID RLPKERSCTS RNLEFKLRRS SSTELLNSLD TESKESKPEH RKSLIFRLNY SASDSNDRSR VKSRLKKFIS RRPSLKTLE KGLIKDQIFG SHLHLVCEHE NSTVPQFVRQ CIKAVERRGL EVDGIYRVSG NLATIQKLRV VVNQEEKLNLD DDSQWEDIHV VTGALKMFFR ELPEPLFPYC FFEQFVEAIK IQDNATRIKA VKTLVKKLPR PNYDTMKVLF EHLKKIAAKE SVNLMSTQSL GIVFGPTLLR PEKETGNMAV HMLYQIQIVE LMLSEYSKIF GSEED
Source	E.coli
Target Names	ARHGAP15
Protein Names	Recommended name: Rho GTPase-activating protein 15 Alternative name(s): ArhGAP15 Rho-type GTPase-activating protein 15
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.