



# Recombinant Mannonate dehydratase (uxuA)

<b>Product Code</b>	CSB-BP723260TNK
<b>Abbreviation</b>	uxuA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q60040
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Thermotoga neapolitana
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKLVFRWFGE KHDTVTLEQI RQIPGVEGVV GALFDIPVGE VWPLEEIMKL KETVERAGLK LEVIESVNVH EDIKLGLPTR DRYIENYKET IRNLAKAGVK VVCYNFMPVF DWMRTDLHKK LPDGSETMEY DHS LiegVTP DELIERVKEG SEGFVLPWE WDRLEKLRET FELYKNVDEE KLFENLVYFL ERVIPVCEEC DVKLAIHPDD PPWSIFGLPR IITNKENIER LLKAVDSPYN GITLCMGSLG ANPENNIPEM IRYFGKMGRI HFAHVRNLKF TGEKSFYETA HPSFCGSHDL FEVMKAFHDI GYEGYIRPDH GRLIWGEKAR PGYGLYDRAL GATYILGLWE AISKMKERYC
<b>Source</b>	Baculovirus
<b>Target Names</b>	uxuA
<b>Protein Names</b>	Recommended name: Mannonate dehydratase EC= 4.2.1.8 Alternative name(s): D-mannonate hydrolase
<b>Expression Region</b>	1-360
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.