



# Recombinant Mesocricetus auratus Ribonucleoside-diphosphate reductase subunit M2 (RRM2)

<b>Product Code</b>	CSB-YP723320MRG
<b>Abbreviation</b>	RRM2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q60561
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mesocricetus auratus (Golden hamster)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFSVRVPLAT ITDQQQLQVS PLKALSLADK ENTPPSLSAT PVLASKVARR ILQDVAEPES KVSTNPSVED EPLLRENPRR FVVFPIEYHD IWKMYKKA EASFWTAAEEVDL SKDIQHWEAL KPDERHFISH VLAFFAASDG IVNENLVERF SQEVQVTEAR CFYGFQIAME NIHSEMYSL IDTYIKDSKE REYLFNAIET MPCVKKKADW ALRWIGDKEA TYGERVVAFA AVEGIFFSGS FASIFWLKKR GLMPGLTFSN ELISRDEGLH CDFACLMFKH LVHKPSEQRV QEITNAVRI EQEFLTEALP VKLIGMNCTL MKQYIEFVAD RLMLELGFNK IFKVENPFDF MENISLEGKT NFFEKRVGGEY QRMGVMSNSF TLDADF
<b>Source</b>	Yeast
<b>Target Names</b>	RRM2
<b>Protein Names</b>	Recommended name: Ribonucleoside-diphosphate reductase subunit M2 EC=1.17.4.1 Alternative name(s): Ribonucleotide reductase small chain Ribonucleotide reductase small subunit
<b>Expression Region</b>	1-386
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.