



Recombinant Mouse Ski-like protein (Skil)

Product Code	CSB-EP723333MO-B
Abbreviation	Skil
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60665
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MENLQSKFSL VQGSNKKLNG MEDDGSPPVK KMMTDIHANG KTLTKVKKEH LDDYGDASVE PDGEHAKRNR VSLPETLNLN PSLKHTLAQF HLSSQSSLGG PAAFSARYSQ ESMSPVFLP LPSPQVLP GP LLIPSDSSTE LTQTLLEGES ISCFQVGGEK RLCLPQVLNS VLREFSLQQI NTVCEDELYIY CSRCTSDQLH ILKVLGILPF NAPSCGLITL TDAQRLCNAL LRPRTFPQNG SILPAKSSLA QLKETGSAFE VEHECLGKCQ GLFAPQFYVQ PDAPCIQCLE CCGMFAPQTF VMHSHRSPDK RTCHWGFESA KWHCYLHVNQ KYLGTPEEK LKIILEEMKE KFSMRNGKRI QSKTDTPSGM ELPSWYPVIK QEGDHVPQTH SFLHPSYYLY MCDKVVAPNV SLTSAASQSK EATKAETSKS TSKQSEK PHE SSQHKT VSY PDVSL EEQEK MDLKTSRELY SCLDSSISNN STSRKKSESA VC SLVRGTSK RDESDSSPLL VRDGEDDKGK IMEDVMRTYV RQQEKLNSIL QRKQQLQMEV EMLSSSKAMK ELTEEQQNLQ KELESLQSEH AQRMEEFYIE QRDLEKKLEQ VMQQKCTCDS TLEKDREA EY AAQLAELRQR LDHAEADRQE LQDEL RQERE ARQKLEMMIK ELKLQIGKSS KPSKD
Source	E.coli
Target Names	Skil
Protein Names	Recommended name: Ski-like protein Alternative name(s): Ski-related oncogene Ski-related protein
Expression Region	1-675
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.