



Recombinant Mouse Origin recognition complex subunit 2 (Orc2)

Product Code	CSB-MP723350MO
Abbreviation	Orc2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60862
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSTLQLKETK VPSVQFVGDD DVLSHILDRE GGTKLKKEKA QLLVNPQKVI KKADCELEKS DLEVLEDQNY VKVLGRNIQE SLGNLSAKDG RNKVVSFQQR KHPEEMTKLA LELAKTSGKK DPLDSNDPEI TKNIAQKSKG HSTSEKAPLV NNNKTEFLST QPHNLRKRII ASRSHYDSES EYSASSSEDD EEATKDEEED TNVARLSQKS QGQNRLLPAP VSKETLPKKK KRDKASDLVE EYFEAHSSSK VLTSDRTLQR LRRARVDQKT LHNLLRKFVP SFSAEIERLN QQHEKLFHKW MLQLHLGFNI VLYGLGSKRD LLEKFRTTML QDSIHVVING YFPGVSVKSI LNSITEDVLS HVGTFQSVLD QRDWIINRFK EDSSLELFL IHNLDSQMLR GDNSQQILGQ LSSLHNVYLI ASIDHLNAPL MWDHAKQSLY NWLWYETTTY SPYTEETSYE NSLLVKQSGS LPLSSLIHVL RSLTPNARGI FRLLMKFQLD NQDSPSYIGL SFQDFYQQCR EAFLVNSDLT LRAQLTEFRD HKLIRTKKGT DGVEYLLIPV DSGILADFLE KEEEEE
Source	Mammalian cell
Target Names	Orc2
Protein Names	Recommended name: Origin recognition complex subunit 2
Expression Region	1-576
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.