



Recombinant Mouse Krueppel-like factor 3 (Klf3)

Product Code	CSB-EP723361MO
Abbreviation	Klf3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q60980
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MLMFDPVPVK QEAMDPVSVS FPSNYIESMK PNKYGVIIYST PLPDKFFQTP EGLTHGIQVE PVDLTVNKRK SPPAAGGSPS SLKFPSHRRA SPGLSMPSSS PPIKKYSPPS PGVQPFQVPL SMPVMAAAL SRHGIRSPGI LPVIQPVVVQ PVPFMYTSHL QQPLMVSLSE EMDNSNSGMP VPVIESYEKP LLQKKIKIEP GIEPQRTDYY PEEMSPPLMN PVSPQALLQ ENHPSVIVQP GKRPLPVESP DTQRKRRIHR CDYDGCNKVY TKSSHLKAHR RTHHTGEKPYK CTWEGCTWKF ARSDDELTRHF RKHTGIKPFQ CPDCDRSFSR SDHLALHRKR HMLV
Source	E.coli
Target Names	Klf3
Protein Names	Recommended name: Krueppel-like factor 3 Alternative name(s): Basic krueppel-like factor CACCC-box-binding protein BKLF TEF-2
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>