



Recombinant Methylococcus capsulatus

Nitrogenase-stabilizing/protective protein nifW

(nifW)

Product Code	CSB-MP723375MFM
Abbreviation	nifW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60C56
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MNPVLQRLQS FSAAEFLDF FGVEYEPSVV HVNRLHILKR FNQYLNRSVP PDDMDEV TAM ATCKALLKQA HDDFVKSTAA QEKVFKVFQD QDGKSISLDS LKASLATRGQ RA
Source	Mammalian cell
Target Names	nifW
Protein Names	Recommended name: Nitrogenase-stabilizing/protective protein nifW
Expression Region	1-112
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.