



Recombinant *Oryza sativa* subsp. japonica PHD finger protein ALFIN-LIKE 5 (Os05g0419100, LOC_Os05g34640)

Product Code	CSB-EP723380OFG-B
Abbreviation	Os05g0419100, LOC_Os05g34640
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60DW3
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDGGSGGPYT SRTAEEVFRD FRGRRAGMIK ALTTDVEKFY QLCDPEKENL CLYGYNETW EVTLPAEEVP PEIPEPALGI NFARDGMNEK DWLALVAVHS DSWLLAVAFY FAARFGFDKE ARRLFNMIN NLPTIFEVVT GAAKKQTKEK APNSTNKPNK PSSKMQRPE SHSKAPKPPA PPKDDDESGD EYADEEEEEER DNTLCGSCGT NDGKDEFWIC CDSCERWYHG KCVKITPARA EHIKHYKCPD CGNKRARA
Source	E.coli
Target Names	Os05g0419100
Protein Names	Recommended name: PHD finger protein ALFIN-LIKE 5
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.