



Recombinant *Macaca fascicularis* Aspartoacylase (ASPA)

Product Code	CSB-MP723388MOV
Abbreviation	ASPA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60HH2
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	≥85% (SDS-PAGE)
Sequence	MTSCHIAEEP IKKVAIFGGT HGNELTGVFL VKHWLENGAE IQRTGLEVKP FITNPRAVKK CTRYIDCDLN RIFDLENLGK KMSEDLPYEV RRAQEIHFLF GPKDGEDSYD IIFDLHNTTS NMGCTLIED SRNNFLIQMF HYIKTSLAPL PCYVYLIEHP SLKYATTRSI AKCPVGIEVG PQPQGVLRAD ILQMRKMIK HALDFIHFN EGKEFPFCAI EVYKIEKVD YPRDENGEIA AVIHPNLQDQ DWKPLHPGDP MFLTLDGKTM PLGGDCTVYP VVFVNEAAYYE KKEAFAKTTK LMLNAKSIRC SVH
Source	Mammalian cell
Target Names	ASPA
Protein Names	Recommended name: Aspartoacylase EC= 3.5.1.15 Alternative name(s): Aminoacylase-2 Short name= ACY-2
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.