



# Recombinant Mouse 5'-nucleotidase (Nt5e)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP723415MO-B  |
| <b>Abbreviation</b>      | Nt5e  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q61503  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Mus musculus (Mouse)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | WE LTILHTNDVH SRLEQTSDDS TKCLNASLCV GGVARLFTKV<br>QQIRKEEPNV LFLDAGDQYQ GTIWFTVYKG LEVAHFMMNIL GYDAMALGNH<br>EFDNGVEGLI DPLLRNVKFP ILSANIKARG PLAHQISGLF LPSKVLSVGG<br>EVLGIVGYTS KETPFLSNPG TNLVFEDEIS ALQPEVDKLG TLNVNKIIL<br>GHSGFEMDKL IAQKVRGVDI VVGGHSTNFL YTGNNPSKEV PAGKYPFIVT<br>ADDGRQVPVV QAYAFGKYLK YLKVFEFDDKG NVITSYGNPI LLNSSIPEDA<br>TIKADINQWR IKLDNYSTQE LGRTIVYLDG STQTCRFREC NMGNLICDAM<br>INNNLRHPDE MFWNHVSMCI VNGGGIRSPI DEKNNGTITW ENLAAVLPFG<br>GTFDLVQLKG STLKKA FEHS VHRYGQSTGE FLQVGGIHVV YDINRKPWNR<br>VVQLEVLCTK CRVPIYEPL MDKVYKVTLP SYLANGGDGF QMIKDELLKH<br>DSGDQDISVV SEYISKMKVV YPAVEGRIKF S |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | Nt5e  |
| <b>Protein Names</b>     | Recommended name: 5'-nucleotidase Short name= 5'-NT EC= 3.1.3.5<br>Alternative name(s): Ecto-5'-nucleotidase CD_antigen= CD73   |
| <b>Expression Region</b> | 29-551  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |