



# Recombinant Mouse Keratin, type I cuticular Ha3-II (Krt33b)

<b>Product Code</b>	CSB-MP723431MO
<b>Abbreviation</b>	Krt33b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q61897
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPYNCCLPAM SCRTSCSSRP CVPPSCHGCT LPGACNIPAN VGNCNWFCEG SFNGNERETM QFLNDRLAS Y MEKVRQLERE NAELECKIQE RNQQQDPLVC PAYQAYFRTI EELQQKILCG KSENGRLVVQ IDNAKLAADD FRTKYESELS LRQLVESDLN GLRRILDELT LCKSDLEAQV ESLRQELLCL KQNHEEEVNT LRCQLGDRLN VEVDAAPTVD LNRVLNETRC QYEALVETNR REVEEWFTTQ TEELNKQVVS SSEQLQSCQA EIIELRRTVN ALEIELQAQH SMRNSLENTL AESEARYSSQ LSQVQCLITN VESQLGEIRA DLERQNQEYQ VLLDVKARLE CEINTYRGLL ESEDCKLPCN PCATTNACDK PIGPCVTNPC VTRSRCGPCN SFGC
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Krt33b
<b>Protein Names</b>	Recommended name: Keratin, type I cuticular Ha3-II Alternative name(s): Hair keratin, type I Ha3 Keratin, type I cuticular Ha3 Keratin-33B Short name= K33B
<b>Expression Region</b>	1-404
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is a member of the keratin gene family. It is one of the type I hair keratin genes which are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. As a type I hair keratin, it is an acidic protein which heterodimerizes with type II keratins to form hair and nails. There are two isoforms of this protein, encoded by two separate genes, KRTHA3A and KRTHA3B.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.