



# Recombinant Mouse Follistatin-related protein 1 (Fstl1)

<b>Product Code</b>	CSB-YP723459MO
<b>Abbreviation</b>	Fstl1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q62356
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	EE EPRSKSKICA NVFCGAGREC AVTEKGEPTC LCIEQCKPHK RPVCGSNGKT YLNHCELHRD ACLTGSKIQV DYDGHCKEKK SASPSASPVV CYQANRDELR RRLIQWLEAE IIPDGWFSKG SNYSEILDKY FKSFDNGDSH LDSSEFLKFV EQNETAINIT TYADQENNKL LRSLCVDALI ELSDENADWK LSFQEFLKCL NPSFNPPEKK CALEDETYAD GAETEVD CNR CVCSCGHWVC TAMTCDGKNQ KGVQTHTEEE KTG YVQELQK HQGTAEKTKK VNTKEI
<b>Source</b>	Yeast
<b>Target Names</b>	Fstl1
<b>Protein Names</b>	Recommended name: Follistatin-related protein 1 Alternative name(s): Follistatin-like protein 1 TGF-beta-inducible protein TSC-36
<b>Expression Region</b>	19-306
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a protein with similarity to follistatin, an activin-binding protein. It contains an FS module, a follistatin-like sequence containing 10 conserved cysteine residues. This gene product is thought to be an autoantigen associated with rheumatoid arthritis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.