



Recombinant Rat Rab-3A-interacting protein (Rab3ip)

Product Code	CSB-MP723472RA
Abbreviation	Rab3ip
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62739
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MANDPLEGFH EVNLASPTSP DLLGVCDPGT QEQTTSPSVI YRPHPSTLCS ATIQANALNL SDLPTQPVYS SPRHLNCAEI SNISIHVPEP ASSVASEVAA GLTRFTSRKD SCNAEREFLLQ GATVTEASAG NDDIFGLSTD SLSRLRSPSV LEVREKGYER LKEELAKAQR ELKLDDEECE RLSKVRDQLG QELEELTASL FEEAHKMOVRE ANVKQATAEK QLKEAQQKID VLQAEVAALK TLVLSSTSPTS PTQEPLADGK TPFKGGHTRN KSTSSAMSGS HQDFSAIQAI VKDCREADLS LYNEFRSWKD EPTMDRTPCPF LDKIYQEDIF PCLTFKSEL ASAVLEAVEN NTLSIEPVGL QPIRFVKASA VECGGPKKCA LTGQSKPCKH RIKLGDSSSY YYISPVCRYR ITSVCNFFTY IRYIQQGLVK QQDVDQMFWE VMQLRKEMSL AKLGYFKEEL
Source	Mammalian cell
Target Names	Rab3ip
Protein Names	Recommended name: Rab-3A-interacting protein Short name= Rab3A-interacting protein Alternative name(s): Rabin-3 SSX2-interacting protein
Expression Region	1-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.