



Recombinant Bacillus cereus Nucleoid occlusion protein (noc)

Product Code	CSB-BP723551BAAD
Abbreviation	noc
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q630C1
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MKNTFSRLFG FGDKESEFEL QDESHEEIDK KVYEEIQEIP IVNITPNRYQ PRTVFDDARI EELALTIRTH GLIQPIVVRQ YEDDKYEIIA GERRFRAATK LGWEKVPAIL KNLNDTETAS VALIENLQRE ELTAIEEAVA YQKLIELHNL TQEALAQLRG KGQSTIANKL RLLKLPEEIK NALLEKSITE RHARALIPLK NEELQLKVLQ EIVEKQLNVK QTEERITKLL EEAKPKRKAK QKAIVSRDTRI AMNTIRQSLQ MVADSGNLVN SEEEEFDEYY QITIQIPKKK
Source	Baculovirus
Target Names	noc
Protein Names	Recommended name: Nucleoid occlusion protein Short name= Noc
Expression Region	1-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.