



Recombinant Bacillus cereus Glucose-6-phosphate isomerase (pgi)

Product Code	CSB-YP723569BAAD
Abbreviation	pgi
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q632G4
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MSTHVTFDYS KALSFIDGEHE ITYLRDAVKV THHAIHEKTG AGNDFLGWVD LPLQYDKKEEF ARIQKCAEKI KNDSDILLVV GIGGSYLGAR AAIEMLNHSF YNTLSKEQRK TPQVLFVGQN ISSTYMKDLM DVLEGKDFSI NVISKSGTTT EPALAFRIFR KLLEEKYGKE EARKRIYATT DKARGALKTL ADNEGYESFV IPDDVGGRFS VLTPVGLLPI AVSGLNIEEM MKGAAAGRDD FGTSELEENP AYQYAVVRNA LYNKGKTIEM LVNYEPALQY FAEWWKQLFG ESEGKDQKGI FPSSANFSTD LHSLGQYVQE GRRDLFETVL KVGKSTHELT IESEENDLDG LNYLAGETVD FVNTKAYEGT LLAHSDGGVP NLIVNIPELN EYTFGYLVYF FEKACAMSGY LLGVNPFDPQ GVEAYKKNMF ALLGKPGFEE LKAELEERLK
Source	Yeast
Target Names	pgi
Protein Names	Recommended name: Glucose-6-phosphate isomerase Short name= GPI EC= 5.3.1.9 Alternative name(s): Phosphoglucose isomerase Short name= PGI Phosphohexose isomerase Short name= PHI
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.