



Recombinant Bacillus cereus Ribonuclease Z (rnz)

Product Code	CSB-BP723591BAAD
Abbreviation	rnz
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q635E2
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MEFVFLGTGA GVPSKGRNVS AIALQLLEER GQTWLFDCGE ATQHQILHTS VRPRRIEKIF ITHLHGDHIF GLPGLLGSR S FQGGTTPLTV YGPKGIIKQFI EVALSVSTTH VKYPLEIVEI TEEGTVFEDN EFHVETKRLS HGIECFGYRI IEKDIQGALL VDKLLEMGVK PGPLFKRLKD GEVVELENGT ILNGKDFIGP PQKGRVITIL GDTRYCEASR ELAQDADLVV HEATFAAEDE QQAYDYFHST SKQAASIALQ ANAKRLILTH ISSRYQGDTY KELLKEAREL FSNTEIATDL KSFPVDR
Source	Baculovirus
Target Names	rnz
Protein Names	Recommended name: Ribonuclease Z Short name= RNase Z EC= 3.1.26.11 Alternative name(s): tRNA 3 endonuclease tRNase Z
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>