



Recombinant Mouse RNA-binding protein Raly (Raly)

Product Code	CSB-YP723712MO
Abbreviation	Raly
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64012
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SLKIQTSNV TNKNDPKSIN SRVFIGNLNT AVVKKSDVET IFSKYGRVAG CSVHKGAFV QYANERHARA AVLGENGRVL AGQTLINMA GEPKPNRPKG LKRAATAIYS GYSFDYDYYQ DYFCARLFDY RGRLSPVPVP RAVPVKRPRV TVPLVRRVKT TIPVKLFARS TAVTTGSAKI KLKSSSELQTI KTELQIKSN IDALLGRLEQ IAEEQKANPD GKKKGDSSSG GGGSSSGGGG SSNVGGGSSG GSGSCSSSSR LPAPQEDTAS EAGTPQGEVQ TRDDGDEEGL LTHSEEELEH SQDTDAEDGA LQ
Source	Yeast
Target Names	Raly
Protein Names	Recommended name: RNA-binding protein Raly Alternative name(s): Maternally-expressed hnRNP C-related protein hnRNP associated with lethal yellow protein
Expression Region	2-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	In infectious mononucleosis, anti-EBNA-1 antibodies are produced which cross-react with multiple normal human proteins. The cross-reactivity is due to anti-gly/ala antibodies that cross-react with host proteins containing configurations like those in the EBNA-1 repeat. One such antigen is RALY which is a member of the heterogeneous nuclear ribonucleoprotein gene family.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.