



Recombinant Mouse Basic helix-loop-helix transcription factor scleraxis (Scx)

Product Code	CSB-MP723719MO
Abbreviation	Scx
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64124
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MSFAMLR SAP PPGRYLYPEV SPLSEDED RG SESSGSDEKP CRVHAARCGL QGARRRAGGR RAAGSGPGPG GRPGREPRQR HTANARERDR TNSVNTAFTA LRTLIPTEPA DRKLSKIETL RLASSYISHL GNVLLVGEAC GDGQPCHSGP AFFHSGRAGS PLPPPPPPPP LARDGGENTQ PKQICTFCLS NQRKLSKDRD RKTAIRS
Source	Mammalian cell
Target Names	Scx
Protein Names	Recommended name: Basic helix-loop-helix transcription factor scleraxis
Expression Region	1-207
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.