



Recombinant Guinea pig Properdin (CFP)

Product Code	CSB-BP723721GU
Abbreviation	CFP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64181
Product Type	Recombinant Protein
Immunogen Species	Cavia porcellus (Guinea pig)
Purity	≥85% (SDS-PAGE)
Sequence	DPVL CFSQYEDSSS KCKDLLGKDV SLEDCCLNAA YAFKERDNGH CQACRSPRWS PWSSWAPCSV SCSEGSQLRH RRCIGWGGQC SENKAPGTLE WQLQACEEQQ CCPEMGGWSN WEPWGPCTVT CSKGTRIRRR VCNNPAPKCG GHCPGVAQES EACDTQQVCP THGAWGPWGP WSSCLSSCHG GPHKPVETRS RTCSAPEPSK NPPGNPCPGT AYEQQSCAGL PPCPVAGGWG PWGSVSPCSV TCGLGQILEQ RKCDNPVPQH GGSFCTGDDT RAHICNTAVP CPVDGEWEPW GDWSTCTRPH LSAIRCKEIV GQQTRVRICK GRKFNGQRCP GKHQEIRHCY NIQNCIFGEK GSWSQWTPWG LCTPPCGANP TRVRQRRCMA SLPKFSPTVS VVEGQGEKNV TFWGKPLAQC EELQGQKVLL EEKRPCLHVP ACKDPEEEEEP
Source	Baculovirus
Target Names	CFP
Protein Names	Recommended name: Properdin Alternative name(s): Complement factor P
Expression Region	27-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.