



Recombinant Rat DALR anticodon-binding domain-containing protein 3 (Dalrd3)

Product Code	CSB-MP723728RA
Abbreviation	Dalrd3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q641Y9
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MATGRLGVGE TLQALNAAVG PGSPVWFKET HARHLRVRDF LAPRSALQAR FRDGGVPECV FRAVSCCLQGP GVAPVLRCAP TPAGLSLQLQ RPAVFEHVLG ALASYATPAK PASPGPRVIL HCPALRCNPD TLRLSQLRAV LVADHLVRVL RAHGVRVCTV PPVRDPHMST FLQKLREVEWP SASKNTSTET LRTCVLAKFN VSKEETLPPG VLGRCLCKEL VEQRHAAGYD PSIDHCLVTE DLLSVLSELQ EAARHWPEDG HPAPVGGPDA GVDDCVVIHV VSCEEAFQQQ KLDLLWQKVD DQVPHRQKHL VCGPVKVAGV PGTQMTAPEY HRLRHAQVCK ASAQKHGGDL AQDPAWTETF DILSVATIKF EMLSTAPQQG LLAHSTIST KGTKSGTFVM YNCARLATLF EGYKHGMEQG LYPTFPLVSS LDFSLLHDEG EWLLLFNSVL PFLDLLSQTV SLASTPGLHI PVRTEMVCKF LVQLSMDFSS YYNRVHILGE PRPHLFGQMF ARLQLLRAVR EVFHTGLAML GLPPLSHI
Source	Mammalian cell
Target Names	Dalrd3
Protein Names	Recommended name: DALR anticodon-binding domain-containing protein 3
Expression Region	1-538
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.