



# Recombinant Rat Heparin cofactor 2 (Serpind1)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP723731RA   |
| <b>Abbreviation</b>      | Serpind1   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q64268   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Rattus norvegicus (Rat)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | AEQLTKE NLTVSLLPPN FHKENTVTND WIPEGEEDDD YLDLEKLLSE<br>DDDIYVVDVA VSPTDSESSA GNILQLFQ GK SRIQRLNILN AKFAFNLYRV<br>LKDQATSSDN IFIAPVGIST AMGMISLGLR GETHEEVHSV LHKDFVNAS<br>SKYEVTTIHN LFRKLTHRLF RRNFGYTLQS VNDLYIQKQF PIREDFKAAM<br>REFYFAEAQE ADFSDPAFIS KANSHILKLT KGLIKEALEN TDSATQMMIL<br>NCIYFKGAWM NKFPVEMTHN HNFRLNEREV VKVSMMQTKG NFLAANDQEL<br>DCDILQLEYV GGISMLIVIP RKLSGMKTLE AQLTPQVVER WQKSMTNRTR<br>EVLLPKFKLE KNYNLVEVLK SMGITKLFNK NGNMSGISDQ RIIIDLFKHQ<br>STITVNEEGT QAAAVTTVGF MPLSTQVRFT VDRPFLFLVY EHRTSCLLFM<br>GRVANPAKS |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | Serpind1   |
| <b>Protein Names</b>     | Recommended name: Heparin cofactor 2 Alternative name(s): Heparin cofactor II Short name= HC-II Protease inhibitor leuserpin-2 Serpin D1   |
| <b>Expression Region</b> | 24-479   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length of Mature Protein  |
| <b>Target Details</b>    | The product encoded by this gene is a serine proteinase inhibitor which rapidly inhibits thrombin in the presence of dermatan sulfate or heparin. The gene contains five exons and four introns. This protein shares homology with antithrombin III and other members of the alpha 1-antitrypsin superfamily. Mutations in this gene are associated with heparin cofactor II deficiency.   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final   |



concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.