



Recombinant Mouse Sequestosome-1 (Sqstm1)

Product Code	CSB-MP723738MO
Abbreviation	Sqstm1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64337
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	ASFTVKAYL LGKEEATREI RRFSCFCFSPE PEAEAQAAAG PGPCERLLSR VAVLFPTLRP GGFQAHYRDE DGDLVAFSSD EELTMAMSYV KDDIFRIYIK EKKECRREHR PPCAQEAPRN MVHPNVICDG CNGPVVGTRY KCSVCPDYDL CSVCEGKGLH REHSLKIFPN PFGHLSDSFS HSRWLRKCLKH GHFGWPGWEM GPPGNWSPRP PRAGDGRPCP TAESASAPPE DPNVNFKNV GESVAAALSP LGIEVDIDVE HGGKRSRLTP TPESSTGT EDKSNTQPSS CSSEVSKPDG AGEQPAQSLT EQMKKIALES VGQPEEQMES GNCSGGDDDW THLSSKEVDP STGELQSLQM PESEGPSSLD PSQEGPTGLK EAALYPHLP EADPRLIESL SQMLSMGFSD EGGWLTRLLQ TKNYDIGAAL DTIQYSKHPP PL
Source	Mammalian cell
Target Names	Sqstm1
Protein Names	Recommended name: Sequestosome-1 Alternative name(s): STONE14 Ubiquitin-binding protein p62
Expression Region	2-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a multifunctional protein that binds ubiquitin and regulates activation of the nuclear factor kappa-B (NF-κB) signaling pathway. The protein functions as a scaffolding/adaptor protein in concert with TNF receptor-associated factor 6 to mediate activation of NF-κB in response to upstream signals. Alternatively spliced transcript variants encoding either the same or different isoforms have been identified for this gene. Mutations in this gene result in sporadic and familial Paget disease of bone.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.