



Recombinant Mouse POU domain class 2-associating factor 1 (Pou2af1)

Product Code	CSB-MP723755MO
Abbreviation	Pou2af1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64693
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MLWQKSTAPE QAPAPRPYQ GVRVKEPVKE LLRRKRGHTS VGAAGPPTAV VLPHQPLATY STVGPSCCLDM EVSASTVTEE GTLCAGWLSQ PAPTALQPLA PWTPYTEYVS HEAVSCPYST DMYVQPVCPS YTVVGPSSVL TYASPPLITN VTPRSTATPA VGPQLEGPEH QAPLTYFPWP QPLSTLPTSS LQYQPPAPTL SGPQFVQLPI SIPEPVLQDM DDPRAISSL TIDKLLLEEE ESNTYELNHT LSVEGF
Source	Mammalian cell
Target Names	Pou2af1
Protein Names	Recommended name: POU domain class 2-associating factor 1 Alternative name(s): B-cell-specific coactivator OBF-1 BOB-1 Short name= BOB1 OCA-B OCT-binding factor 1
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.