



Recombinant Bacteroides fragilis Elongation factor Ts (tsf)

Product Code	CSB-MP723786BDQ
Abbreviation	tsf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64P30
Product Type	Recombinant Protein
Immunogen Species	Bacteroides fragilis (strain YCH46)
Purity	≥85% (SDS-PAGE)
Sequence	MAVTMADITK LRKMTGAGMM DCKNALTDAE GDFDKAMKII REKGQAVAAK RSDREASEGC VLVKVEEGFG AIALKCETD FVAQNADSVK LTQDILDAAV ANKCKTLEEV LALPMGDATV AQAVTDRTGI TGEKMELDGY MVLEGATIAA YNHMMNRNGLC TMVAFNKKVD EQLAKQVAMQ VAAMNPIAVD EDGVSEEVKQ KEIEVAVEKT KVEQVQKAVE AALKKANINP AHVDSSEHME SNMAKGWITA EDVAKAKEII ATVSAEKAAN MPEQMIQNIA KGR LAKFLKE VCLLNQEDIM DAKKTVREVL KEADPELKVV DFKRFTLRAE
Source	Mammalian cell
Target Names	tsf
Protein Names	Recommended name: Elongation factor Ts Short name= EF-Ts
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.