



Recombinant Bacteroides fragilis Riboflavin biosynthesis protein ribBA (ribBA)

Product Code	CSB-MP723809BDQ
Abbreviation	ribBA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64YT3
Product Type	Recombinant Protein
Immunogen Species	Bacteroides fragilis (strain YCH46)
Purity	≥85% (SDS-PAGE)
Sequence	MEPIRLNTIE EAIADFKEGN FVIVVDDedr ENegDFIIAA EKITPEKVNf MLTHGRGVLC APITEERCAE LELDMQVSSN TSIYETPFTV TVDLLEGCTT GVSMHDRAMT IRALADPKTK PADLGRPGHI NPLRARSRGV LRRAGHTEAS VDLAKLAGLY PAAALIEIIN EDGTMARLPQ LVEVARRFGL KIISIKDLIA YRLQMESIVD RGVEVDMPTQ FGHFRLIPFR QKSNGMEHIA LIKGTWDTDE PILVRVHSSC MTGDIFGSCR CECGEQLHKA MEMIEAAGKG VIVYMNQEGR GIGLMNKIAA YKLQEEGYDT VDANLHLGFD ADERDYGvGA QILREIGVKK MKLMTNPNVK RIGLEAYGLE ITENVGIEIK PNPYNERYLK TKKDRMGHTL HFNK
Source	Mammalian cell
Target Names	ribBA
Protein Names	Recommended name: Riboflavin biosynthesis protein ribBA Including the following 2 domains: 3,4-dihydroxy-2-butanone 4-phosphate synthase Short name= DHBP synthase EC= 4.1.99.12 GTP cyclohydrolase-2 EC= 3.5.4.25 Alterna
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.