



Recombinant *Oryza sativa* subsp. *japonica* Protein MEI2-like 6 (ML6)

Product Code	CSB-YP723832OFG
Abbreviation	ML6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q652K6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSCLSRLNAS AAPWEPPVAR AMAVEQYCPP PQSLLPPPPL PPVAVPTTCG CAACLQGC FV PVGVQA AAFPH AAAGCAPP PP VMPVMIVYRV VQPPPPAAHA TRCQITEIED GGGVETAKAV DGDEPQPFIR TVRSTRRRKA AAIRLPKA FR AALLPPPPPP CALGFTATTT SLMIRNIPNK FLKARLMAIL DQHCADENGK CHRRGGGGGR SVVKSEYDFF YVPIDFKTGF NKGYAFV NMT TATAARRLRA FLQDHRWDAA MSGKVCDVVP AAIQGLDAFV AHFSASCFPC RTKEFLPVWF EPPRDGEQQT KAHVVGRLVV RPR
Source	Yeast
Target Names	ML6
Protein Names	Recommended name: Protein MEI2-like 6 Short name= OML6 Alternative name(s): MEI2-like protein 6
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.