



Recombinant Phoca vitulina Lysozyme C (LYZ)

Product Code	CSB-YP723845EUR
Abbreviation	LYZ
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q659U1
Product Type	Recombinant Protein
Immunogen Species	Phoca vitulina (Harbor seal)
Purity	≥85% (SDS-PAGE)
Sequence	KI FERCDLARTL KRLGLAGFKG VSLANWMCLA KWESNYNTKA TNYNPGSRST DYGIFQINSR YWCNDGKTPR AVNACHIPCS DLLKDDITQA VACAKRVVSD PNGIRAWVAW RAHCENQDVS QYVRNCGV
Source	Yeast
Target Names	LYZ
Protein Names	Recommended name: Lysozyme C EC= 3.2.1.17 Alternative name(s): 1,4-beta-N-acetylmuramidase C
Expression Region	19-148
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>