



Recombinant Ribonuclease Z (rnz)

Product Code	CSB-MP723976BUI
Abbreviation	rnz
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q660C0
Product Type	Recombinant Protein
Immunogen Species	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia baviensis)
Purity	≥85% (SDS-PAGE)
Sequence	MGFNINIIGT GGTRPLHNR Y LSSVLIEYNG DNFLFDCGEG TQMSLRKQKI SWQKIKMICI THLHADHITG LLGIVMLMSQ SGETRKEPLI IAGPVGIKNY TKTNIDILKI YKNYEIIYKE IIDKTEKII YEDKTKKIEY TRLKHSIECV GYLFIEKDKP GKFNTEKAE LNIPKGPPIRK TLQDGKEILI DGKIIKPSEI LGESKKGLKV AYITDTGYFK ELIQQIKNFN LVIIESTFKN ELKKEADKKL HLTAGGAANI VKQAKVLKTG LIHFSERYTL RKDLENLLKE AKLEYPDGEI FLTKDGMRL E ANKNNFIIK
Source	Mammalian cell
Target Names	rnz
Protein Names	Recommended name: Ribonuclease Z Short name= RNase Z EC= 3.1.26.11 Alternative name(s): tRNA 3 endonuclease tRNase Z
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>