



# Recombinant GTPase Der (der)

<b>Product Code</b>	CSB-MP723983BU1
<b>Abbreviation</b>	der
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q661B2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBI) (Borrelia baviensis)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>           MLSYKKVLIV GRPNVGKSAL FNRILDAKRS ITESTYGVTR DLVEEVCKVG            SFNFKLIDTG GFTILKDEIS KIVVQKVLSS LEKVDLILLV LDVNEILLED            YEIIERLRKY SSKVILVLNK VDTKDKFLA HKFHNLGFKR YFLVSALHRR            GITKLRDFLK VEVGRVNIEE EVNIKVGIIIG KPNSGKSTLI NYLSGNEISI            VSDQPGTTRD FIKTKLTRNG KVFEIIDTAG IRRRARVNEV VEYYSVNRAL            KVIDIVDIVF LLIDVEEELT SQDKKIAHYV TTKGKGIIIV FSKWDLLEES            KGYFETLKGR VKFFFPVLNF APIFRISVYK KIGLDSLDFKE ALKVKDQLEL            KTNTPDNLKM LNLWIKDYHL NISHKIKYIT QVSTNPVKFI LFANKIKNFP            NSYYNYLVNN LRKIGYKNIP ILVELREKIR DLK         </p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	der
<b>Protein Names</b>	Recommended name: GTPase Der Alternative name(s): GTP-binding protein EngA
<b>Expression Region</b>	1-433
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.