



Recombinant Ribosomal RNA small subunit methyltransferase G (rsmG)

Product Code	CSB-MP723996BUI
Abbreviation	rsmG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q662I7
Product Type	Recombinant Protein
Immunogen Species	Borrelia bavariensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia bavariensis)
Purity	≥85% (SDS-PAGE)
Sequence	MISDIEFAFS ECNFQFAYKD LQKINLYIKR ILLLNTRFNL ISNSNSNFNS ILNLHVIDSL LGLPTIKEIN PSEILDVGSG SGFPGIVLAI FDTSRKYLL ERSKKKSTFL KMIKLELDLE NVKILEYEIE REKKKYEFIT IRAFRSMNEY ALVLKNLLKN GGLIMAYKGGK FDKINLEVNQ IKDLFSKIEV KSLNSKLSLD RNLVLLYR
Source	Mammalian cell
Target Names	rsmG
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase G EC=2.1.1.- Alternative name(s): 16S rRNA 7-methylguanosine methyltransferase Short name= 16S rRNA m7G methyltransferase
Expression Region	1-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.