



Recombinant Yersinia pseudotuberculosis serotype I Maltoporin 2 (IamB2)

Product Code	CSB-EP724028YAH
Abbreviation	IamB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q666Z8
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I
Purity	>85% (SDS-PAGE)
Sequence	IEKIDFH GYMRAVGVS SDGGLAEWQK TMVGRLGNES DTYGEIGLGA EVYKKEDVSF YLESMVSMLS DGSNDSETTI GDDAQFGLRQ LNLQIKGLIP GDKEAVIWGG KRYYQRHDLH IIDTKYWNIS GSGAGIENYT VGPGAVSVAW VRGDANDVDT RITGDSVDNI NYIDVRYAGF KPWAGSWTEV GIDYAMPNPT KQQKEYGGLY DADNAVMLTG EISQDMFGGY NKLVLQYANK GLAQNMISSQ GGWYDMWHKT DEAKGYRVIN TGLIPITDKF SFNHVLTWGS ANDITEYTDK TNLISLVGRA QYQFTQYVRA IGEVGGFYQK DTYHNGSNYK QGGEKYTIAL GLAEGPDFLS RPELRVFASY LNDSSENGKPF EDGTSNDTWN FGVQVEAWW
Source	E.coli
Target Names	IamB2
Protein Names	Recommended name: Maltoporin 2 Alternative name(s): Maltose-inducible porin 2
Expression Region	24-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.