



Recombinant ProP effector (proQ)

Product Code	CSB-MP724048YAH
Abbreviation	proQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q669V5
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MENQPKLNSS KEVIAFLAER FPLCFTAEGE ARPLKIGIFQ DLVERVQGEE NLSKTQLRSA LRLYTSSWRY LYGVKVGAE RVDLDGNPCGV LEEQHVEHAR KQLEEAKARV QAQRAEQQAK KREAAIAAGE TPEPRRPRPA GKKPAPRREA GAAPENRKPR QSPRPQQVRP PRPQVEENQP RPVPVTDISK LQIGQEIKVR AGKSAMDATV LEIAKDGV RV QLSSGLAMIV RAEHLQF
Source	Mammalian cell
Target Names	proQ
Protein Names	Recommended name: ProP effector
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.