



Recombinant Chaperone surA (surA)

Product Code	CSB-YP724087YAH
Abbreviation	surA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66EQ7
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	APQEVDKVAA VVDNGVVLQS DIDGLLQSVK MNAQQSGQQV PDDSTLRHQI LERLIMDNIQ LQMAKKMGIT ITDQALDKAI ADIAAQNRMT LAQMRSRLAA DGLSYDTYPE QIRKEMLTSE VRNNEVRRRI TILPQEVESL AKQMGNQVSG DTELNLSHIL IPLPENPTQQ QVDQAEDLAN KLVADIKGGA DFGKLAIAN ADSQALKGGQ MGWGKLQELP SLFAERLQSA HKGEIVGPIR SGVGFHILKV NDMRGADQTI SVTEVNARHI LLKPSPMMTD EQARAKLEAA AAEIKSGKTS FATIAKEISQ DPGSAMQGGG LGWASPDIYD PAFRDALMKL KKGEISAPVH SSFGWHLIQL VDTRQVDKTD AAQKERAYRM LFNRFKFAEEA QTWMQEQRAA AYYVKILDGSN AQPQ
Source	Yeast
Target Names	surA
Protein Names	Recommended name: Chaperone surA Alternative name(s): Peptidyl-prolyl cis-trans isomerase surA Short name= PPlase surA EC= 5.2.1.8 Rotamase surA
Expression Region	21-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.