



Recombinant Phosphopentomutase (deoB)

Product Code	CSB-BP724089YAH
Abbreviation	deoB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66EV8
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MKRTFIMVLD SFGIGASADA KKFGDEGADT LGHIAEACAR GEANVGRSGP LTLPNLSRLG LGKAAEESTG TFPVGLDKNA DIIGAYGYAS ELSSGKDTPS GHWEIAGVPV LFDWGYFSDV ENSFPQELLD KLVKRANLPG YLGNCHSSGT VILDQLGEEH MKTGKPIFYT SADSVFQIAC HEETFGLDRL YELCEIAREE LTDGGYNIGR VIARPFIGDK PGHFQRTGNR HDLAVEPPAP TMLKKLVDEK GGEVVSIGKI ADIYAQVGIT QKVKATGLDA LFDATIEEMK KAGDNTIVFT NFVDFDSSYG HRRDVAGYAA ALELFDRLRP ELMALVKEDD ILILTADHGC DPTWPGTDHT REHIPVLVYG PKVKPGSLGH RETFADIGQT VAAAYFGLSPM DYGKNML
Source	Baculovirus
Target Names	deoB
Protein Names	Recommended name: Phosphopentomutase EC= 5.4.2.7 Alternative name(s): Phosphodeoxyribomutase
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.