



Recombinant L-rhamnose isomerase (rhaA)

Product Code	CSB-MP724096YAH
Abbreviation	rhaA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66FF4
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MTNSIEQAWD LAKQRFAAVG VDVDAALTRL DTLPVSMHCW QGDDVTGFED PDGVLTTGGIQ ATGNYPGKAR NATELRSDLE LALALIPGPK RLNLHAIYLE SDTPVARNKI EPRHFSHWVT WAKKHQLGLD FNPSCFSHPL SADGFTLSHA DPEIRQFWIE HCQASRRVSA YFGEQLGTPS VMNIWIPDGM KDTPIDRLAP RQRLLSALDE VISEKLNPAH HIDAVESKLF GIGAESYTVG SNEFYMGYAA SRQTALCLDA GHFHPTEVIS DKISSAMLYV PRLLHVSRLP VRWDSHDVVL LDDETQAIAS EIIRHNLFDL VHIGLDFFDA SINRIAAWVI GTRNMKKALL RALLEPTDML RQLELRGDYT ARLALLEEQK SLPWQAIWEG YCQRNDVPVD ARWLDAVREY EQQILSQR
Source	Mammalian cell
Target Names	rhaA
Protein Names	Recommended name: L-rhamnose isomerase EC= 5.3.1.14
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.